



Triumphant Labs
4858 W Pico Blvd #1031 Los Angeles, CA 90019
hi@triumphantlabs.com
http://triumphantlabs.com

SAMPLE NAME: GLP-3 10mg
SAMPLE CODE: F2W5-M5RF
RECEIVED DATE: 09/14/2026
PUBLISHED DATE: 09/17/2026



DECLARED IDENTITY: Retatrutide
MATRIX TYPE: PEPTIDE
SAMPLE NAME: GLP-3 10mg
SAMPLE SIZE: 10.0 mg
LOT CODE: YR_GL3-10_2609



Variance
3 Vials




9.37 mg
Total Quantity



99.98%
Total Blend Purity

Sample	Test	Specification	Result	Status
Retatrutide	IDENTITY	Retatrutide	Conforms 3/3	Conforms 3/3
Retatrutide	QUANTITY	MEASURE	mean 9.37 mg · SD 0.18 · %RSD 1.88% · n=3	MEASURED
Retatrutide	PURITY	> 98.0%	mean 99.98% · SD 0.00 · %RSD 0.00% · n=3	CONFORMS

HPLC Chromatogram Report



RESULTS REVIEWED BY:

Anthony Burke PhD
Lab Manager
09/17/2026

RESULTS CERTIFIED BY:

Dr. Levi Fried MD
Lab Director
09/17/2026

ACCUMARKLABS.COM
F2W5-M5RF
09/17/2026





Triumphant Labs
4858 W Pico Blvd #1031 Los Angeles, CA 90019
hi@triumphantlabs.com
http://triumphantlabs.com

SAMPLE NAME: GLP-3 10mg
SAMPLE CODE: F2W5-M5RF
RECEIVED DATE: 09/14/2026
PUBLISHED DATE: 09/17/2026

VIAL	Retatrutide		
	Identity	Quantity mg	Purity %
● Vial 1	Retatrutide	9.357 mg	99.983%
● Vial 4	Retatrutide	9.547 mg	99.98%
● Vial 5	Retatrutide	9.196 mg	99.983%

RESULTS CERTIFIED BY:


Anthony Burke PhD
Lab Manager
09/17/2026

RESULTS CERTIFIED BY:


Dr. Levi Fried MD
Lab Director
09/17/2026

ACCUMARKLABS.COM
F2W5-M5RF
09/17/2026





Triumphant Labs
4858 W Pico Blvd #1031 Los Angeles, CA 90019
hi@triumphantlabs.com
http://triumphantlabs.com

SAMPLE NAME: GLP-3 10mg
SAMPLE CODE: F2W5-M5RF
RECEIVED DATE: 09/14/2026
PUBLISHED DATE: 09/17/2026

Endotoxin USP85 LAL

Test	Specification (EU/mL)	Result (EU/mL)	LOQ (EU/mL)	Verdict
Endotoxin USP85 LAL	≤5	< LOQ	2.5	Conforms

Method: Endotoxin testing performed using Limulus Amebocyte Lysate assay in accordance with USP <85> under validated laboratory conditions.

Rapid Sterility Screening (PCR)

Test	Specification	Result	Unit	Verdict
Rapid Sterility Screening (PCR)	Not Detected	Not Detected		Conforms

Method: Rapid sterility screening performed by real-time polymerase chain reaction (qPCR) targeting conserved microbial DNA sequences under validated laboratory conditions.

RESULTS CERTIFIED BY:


Anthony Burke PhD
Lab Manager
09/17/2026

RESULTS CERTIFIED BY:


Dr. Levi Fried MD
Lab Director
09/17/2026

ACCUMARKLABS.COM
F2W5-M5RF
09/17/2026

