




**Triumphant Labs**  
4858 W Pico Blvd #1031 Los Angeles, CA 90019  
hi@triumphantlabs.com  
http://triumphantlabs.com


SAMPLE NAME: Bacteriostatic Water  
SAMPLE CODE: 73LB-JG8R  
RECEIVED DATE: 09/14/2026  
PUBLISHED DATE: 09/21/2026



|                   |                      |
|-------------------|----------------------|
| DECLARED IDENTITY | Bacteriostatic Water |
| MATRIX TYPE:      | BACTERIOSTATIC WATER |
| SAMPLE NAME:      | Bacteriostatic Water |
| SAMPLE SIZE:      | 10 mL                |
| LOT CODE:         | KEI_BW_2609          |

| Sample                      | Test        | Specification   | Result       | Status   |
|-----------------------------|-------------|-----------------|--------------|----------|
| Benzyl Alcohol Assay (HPLC) | HPLC        | 0.9% (v/v) ±20% | 0.87 % (v/v) | CONFORMS |
| Fill volume / Net content   | Gravimetric | —               | 9.86 mL      | MEASURED |
| pH Determination            | pH          | 4.5 – 7.0       | 5.00         | CONFORMS |

## HPLC Chromatogram Report



RESULTS REVIEWED BY:

Anthony Burke PhD  
Lab Manager  
09/21/2026

RESULTS CERTIFIED BY:

Dr. Levi Fried MD  
Lab Director  
09/21/2026

ACCUMARKLABS.COM  
73LB-JG8R  
09/21/2026





Triumphant Labs  
4858 W Pico Blvd #1031 Los Angeles, CA 90019  
hi@triumphantlabs.com  
http://triumphantlabs.com

SAMPLE NAME: Bacteriostatic Water  
SAMPLE CODE: 73LB-JG8R  
RECEIVED DATE: 09/14/2026  
PUBLISHED DATE: 09/21/2026

## Endotoxin USP85 LAL

| Test                | Specification (EU/mL) | Result (EU/mL) | LOQ (EU/mL) | Verdict  |
|---------------------|-----------------------|----------------|-------------|----------|
| Endotoxin USP85 LAL | ≤0.25                 | < LOQ          | 0.2         | Conforms |

Method: Endotoxin testing performed using Limulus Amebocyte Lysate assay in accordance with USP <85> under validated laboratory conditions.

## Rapid Sterility Screening (PCR)

| Test                            | Specification | Result       | Unit | Verdict  |
|---------------------------------|---------------|--------------|------|----------|
| Rapid Sterility Screening (PCR) | Not Detected  | Not Detected |      | Conforms |

Method: Rapid sterility screening performed by real-time polymerase chain reaction (qPCR) targeting conserved microbial DNA sequences under validated laboratory conditions.

RESULTS CERTIFIED BY:

 Anthony Burke PhD  
Lab Manager  
09/21/2026

RESULTS CERTIFIED BY:

 Dr. Levi Fried MD  
Lab Director  
09/21/2026

ACCUMARKLABS.COM  
73LB-JG8R  
09/21/2026

